

Naloxone is a safe, easy-to-use, life-saving medication that has one job: Reverse an opioid overdose by blocking the effects of opioids. Learn how to use naloxone, visit: bit.ly/naloxone-trainpa.

The Pennsylvania Department of Health's Acting Secretary signed a [standing order prescription](#) for naloxone. This means you do not need to get a prescription for naloxone from a healthcare provider. The following products are available under Pennsylvania's standing order prescription and are approved by the Food and Drug Administration (FDA) for use by the public. Visit pa.gov/opioids to access the standing order and learn more.

INTRA-NASAL NALOXONE

Prepackaged intra-nasal naloxone products require little to no assembly and are sprayed into one nostril while a person lays on their back. For people who have no formal training, nasal spray naloxone may be easier to use.

NARCAN Nasal Spray

NARCAN Nasal Spray is needle free and requires no assembly. NARCAN Nasal Spray comes in a single 4 mg dose. Two individually packaged single-spray devices are contained in one NARCAN Nasal Spray package.

Find more information and safety instructions at bit.ly/NARCANinfo



KLOXXADO Nasal Spray

KLOXXADO Nasal Spray is needle free and requires no assembly. KLOXXADO Nasal Spray comes in a single 8 mg dose. Two individually packaged single-spray devices are contained in one package.

Find more information and safety instructions at bit.ly/KLOXXADOinfo.



Luer-Lock Syringes and Mucosal Atomization Devices (MAD)

Luer-Lock Syringes and MAD Nasal Device packaging includes two 2 mL Luer-Jet luer-lock syringes pre-filled with naloxone and two mucosal atomization devices (nasal caps). The nasal cap connects to the syringe. While the nasal cap is pressed against the nostril, the syringe plunger is pressed to release the naloxone.

Find more information and safety instructions at bit.ly/MADNASALinfo.

INTRA-MUSCULAR NALOXONE

PRE-FILLED SYRINGE NALOXONE

Pre-filled injectable syringe naloxone products are injected with a needle into the outer thigh, similar to an Epi-pen. Syringe naloxone products may be easier for people who have experience administering medicine with a needle.

ZIMHI

ZIMHI is an autoinjector with a syringe. The syringe is prefilled with a single 5 mg dose of naloxone. Two single-use syringes are contained in one ZIMHI package.

Find more information and safety instructions at bit.ly/ZIMHIinfo.



NON-PRE-FILLED SYRINGE WITH NALOXONE VIAL

Intra-muscular naloxone can also consist of a non-pre-filled syringe with a vial that contains liquid naloxone. The syringe is used to draw the liquid naloxone into the syringe before administration. Once filled with the correct amount of liquid naloxone, the syringe needle is inserted and administered into the muscle of the upper arm or thigh. Non-pre-filled needle-syringe naloxone products may be easier for people who have experience using a syringe to draw liquid medicine from a vial and administering medicine with a needle.



The Pennsylvania Standing Order for naloxone allows eligible persons to be provided with the following: 1) Naloxone 0.4mg/mL in 1mL single dose vials (quantity to dispense: two single dose vials), 2) one 3 cc, 23g or 25g, 1-inch or 1 ½-inch syringe per dose and 3) patient instructions.

Please ask your pharmacist or healthcare provider for more information and safety instructions for the non-pre-filled syringe naloxone products provided to you. The Department does not endorse a particular manufacturer for the supply of non-pre-filled syringes or vials.